URIN-1 URIN-2 URIN-3 URIN-5 URIN-10 URIN-11

For rapid detection of multiple analytes in human urine. For in vitro diagnostic use only. INTENDED USE

The Urinalysis Reagent Strips (Urine) are firm plastic strips onto which several separate reagent areas are affixed. The test is for the detection of one or more of the following analytes in urine: Ascorbic acid, Glucose, Bilirubin, Ketone (Acetoacetic acid), Specific Gravity, Blood, pH, Protein, Urobilinogen, Nitrite and Leukocytes.

SUMMARY
Urine undergoes many changes during states of disease or body dysfunction before blood composition is altered to a significant extent. Urinalysis is a useful procedure as an indicator of health or disease, and as such, is a part of routine health screening. The Urinalysis Reagent Strips (Urine) can be used in general evaluation of health, and aids in the diagnosis and monitoring of metabolic or systemic diseases that affect kidney function, endocrine disorders and diseases or disorders of the urinary tract.

PRINCIPLE AND EXPECTED VALUES

Ascorbic acid: This test involves decolorization of Tillmann's reagent. The presence of ascorbic acid causes the color of the test field to change from blue-green to orange. Patients with adequate diet may excrete 2-10 mg/dL daily. After ingesting large amounts of ascorbic acid, levels can be around 200 mg/dL.

Glucose: This test is based on the enzymatic reaction that occurs between glucose oxidase, peroxidase and chromogen. Glucose if first oxidized to produce gluconic acid and hydrogen peroxide in the presence of glucose oxidase. The hydrogen peroxide reacts with potassium iodide chromogen in the presence of peroxidase. The extent to which the chromogen is oxidized determines the color which is produced, ranging from green to brown. Low amounts of glucose are normally excreted in urine.3 Glucose concentrations as low as 100 mg/dL, read at either 10 or 30 seconds, may be considered abnormal if results are consistent. At 10 seconds, results should be interpreted qualitatively. For semi-quantitative results, read at 30 seconds only.

Bilirubin: This test is based on azo-coupling reaction of bilirubin with diazotized dichloroaniline in a strongly acidic medium. Varying bilirubin levels will produce a pinkish-tan color proportional to its concentration in urine. In normal urine, no bilirubin is detectable by even the most sensitive methods. Even trace amounts of bilirubin require further investigation. Atypical results (colors different from the negative or positive color blocks shown on the color chart) may indicate that bilirubin-derived bile pigments are present in the urine specimen, and are possibly masking the bilirubin reaction.

Ketone: This test is based on ketones reacting with nitroprusside and acetoacetic acid to produce a color change ranging from light pink for negative results to a darker pink or purple color for positive results. Ketones are normally not present in urine. Detectable ketone levels may occur in urine during physiological stress conditions such as fasting, pregnancy and frequent strenuous exercise.4-6 In starvation diets, or in other abnormal carbohydrate metabolism situations, ketones appear in the urine in excessively high concentration before serum ketones are elevated.7

Specific Gravity: This test is based on the apparent pKa change of certain pretreated polyelectrolytes in relation to ionic concentration. In the presence of an indicator, colors range from deep blue-green in urine of low ionic concentration to green and yellow-green in urine of increasing ionic concentration. Randomly collected urine may vary in specific gravity from 1.003-1.040. Twenty-four hour urine from healthy adults with normal diets and fluid intake will have a specific gravity of 1.016-1.022.8 In cases of severe renal damage, the specific gravity is fixed at 1.010, the value of the glomerular filtrate.

Blood: This test is based on the peroxidase-like activity of hemoglobin which catalyzes the reaction of cumene-hydroperoxide and 3,3',5,5'-tetramethylbenzidine. The resulting color ranges from orange to dark blue. Any green spots or green color development on the reagent area within 60 seconds is significant and the urine specimen should be examined further. Blood is often, but not invariably, found in the urine of menstruating females.

pH: This test is based on a double indicator system which gives a broad range of colors covering the entire urinary pH range. Colors range from orange to yellow and green to blue. The expected range for normal urine specimens from newborns is pH 5-7. The expected range for other normal urine specimens is pH 4.5-8, with an average result of pH 6.

Protein: This reaction is based on the phenomenon known as the "protein error" of pH indicators where an indicator that is highly buffered will change color in the presence of proteins (anions) as the indicator releases hydrogen ions to the protein. At a constant pH, the development of any green color is due to the presence of protein. Colors range from yellow to yellow-green for negative results and green to green-blue for positive results. 1-14 mg/dL of protein may be excreted by a normal kidney.9 A color matching any block greater than trace indicates significant proteinuria. For urine with high specific gravity, the test area may most closely match the trace color block even though only normal concentrations of protein are present. Clinical judgment is required to evaluate the significance of trace results.

Urobilinogen: This test is based on a modified Ehrlich reaction between p-diethylaminobenzaldehyde and urobilinogen acid in strongly acidic medium to produce a pink color. Urobilinogen is one of the major compounds produced in heme synthesis and is a normal substance in urine. The expected range for normal urine with

this test is 0.2-1.0 mg/dL (3.5-17 mg/L). A result of 2.0 mg/dL (35 mg/L) may be of clinical significance, and the patient specimen should be further evaluated.
Nitrite: This test depends upon the conversion of nitrate to nitrite by the action of Gram negative bacteria in the urine. In an acidic medium, nitrite in the urine reacts with p-arsanilic acid to form a diazonium compound. The diazonium compound in turn couples with 1 N-(1-naphthyl)- ethylenediamine to produce a pink color. Nitrite is not detectable in normal urine.3 The nitrite area will be positive in some cases of infection, depending on how long the urine specimens were retained in the bladder prior to collection. Retrieval of positive cases with the nitrite test ranges from as low as 40% in cases where little bladder incubation occurred, to as high as approximately 80% in cases where bladder incubation took place for at least 4 hours.

Leukocytes: This test reveals the presence of granulocyte esterases. The esterases cleave a derivatized pyrazole amino acid ester to liberate derivatized hydroxy pyrazole. This pyrazole then reacts with a diazonium salt to produce a beige-pink to purple color. Normal urine specimens generally yield negative results. Trace results may be of questionable clinical significance. When trace results occur, it is recommended to retest using a fresh specimen from the same patient. Repeated trace and positive results are of clinical significance.

REAGENTS AND PERFORMANCE CHARACTERISTICS/REAGENTS AND PERFORMANCE CHARACTERISTICS

Based on the dry weight at the time of impregnation, the concentrations given may vary within manufacturing tolerances. The following table below indicates read times and performance characteristics for each parameter.

Reagent	Read Time	Composition	Description
Ascorbic Acid (ASC)	30 seconds	2,6-dichlorophenolindophenol; buffer and non-reactive ingredients	Detects ascorbic acid as low as 5-10 mg/dL (0.28-0.56 mmol/L).
Glucose (GLU)	30 seconds	glucose oxidase; peroxidase; potassium iodide; buffer; non-reactive ingredients	Detects glucose as low as 50-100 mg/dL (2.5-5 mmol/L)..
Bilirubin (BIL)	30 seconds	2, 4-dichloroaniline diazonium salt; buffer and non-reactive ingredients	Detects bilirubin as low as 0.4-0.8 mg/dL (6.8-13.6 µmol/L).
Ketone (KET)	40 seconds	sodium nitroprusside; buffer	Detects acetoacetic acid as low as 2.5-5 mg/dL (0.25-0.5 mmol/L).
Specific Gravity (SG)	45 seconds	bromthymol blue indicator; buffer and non-reactive ingredients; 55% poly (methyl vinyl ether/maleic anhydride); 25% sodium hydroxide	Determines urine specific gravity between 1.000 and 1.030. Results correlate with values obtained by refractive index method within ±0.005.
Blood (BLO)	60 seconds	3,3',5,5'-tetramethylbenzidine (TMB); cumene hydroperoxide; buffer and non-reactive ingredients	Detects free hemoglobin as low as 0.015-0.062 mg/dL or 5-10 Ery/µL in urine specimens with ascorbic acid content of <50 mg/dL.
pH	60 seconds	methyl red sodium salt; bromthymol blue; non-reactive ingredients	Permits the quantitative differentiation of pH values within the range of 5-9.
Protein (PRO)	60 seconds	tetrabromophenol blue; buffer and non-reactive ingredients	Detects albumin as low as 7.5-20 mg/dL (0.075-0.2 g/L).
Urobilinogen (URO)	60 seconds	p-diethylaminobenzaldehyde; buffer and non-reactive ingredients	Detects urobilinogen as low as 0.2-1.0 mg/dL (3.5-17 µmol/L).
Nitrite (NIT)	60 seconds	p-arsanilic acid; non-reactive ingredients	Detects sodium nitrite as low as 0.05-0.1 mg/dL in urine with a low specific gravity and less than 30 mg/dL ascorbic acid.
Leukocytes (LEU)	120 seconds	derivatized pyrrole amino acid ester; diazonium salt; 32% w/w buffer; non-reactive ingredients	Detects leukocytes as low as 10-25 white blood cells Leu/µL in clinical urine.

The performance characteristics of the Urinalysis Reagent Strips (Urine) have been determined in both laboratory and clinical tests. Parameters of importance to the user are sensitivity, specificity, accuracy and precision. Generally, this test has been developed to be specific for the parameters to be measured with the exceptions of the interferences listed. Please refer to the Limitations section in this package insert.

Interpretation of visual results is dependent on several factors: the variability of color perception, the presence or absence of inhibitory factors, and the lighting conditions when the strip is read. Each color block on the chart corresponds to a range of analyte concentrations.

PRECAUTIONS

- For in vitro diagnostic use only. Do not use after the expiration date.
- The strip should remain in the closed canister until use.
- Do not touch the reagent areas of the strip.
- Discard any discolored strips that may have deteriorated.
- All specimens should be considered potentially hazardous and handled in the same manner as an infectious agent.

The used strip should be discarded according to local regulations after testing.

STORAGE AND STABILITY

Store as packaged in the closed canister either at room temperature or refrigerated (2-30°C). Keep out of direct sunlight. The strip is stable through the expiration date printed on the canister label. Do not remove the desiccant. Remove only enough strips for immediate use. Replace cap immediately and tightly. **DO NOT FREEZE.** Do not use beyond the expiration date.

Note: Once the canister has been opened, the remaining strips are stable for up to 3 months. Stability may be reduced in high humidity conditions.

SPECIMEN COLLECTIONS AND PREPARATION

A urine specimen must be collected if possible. Do not centrifuge. The urine testing cannot be done within an hour immediately and let it return to room temperature. Prolonged storage of unpreserved urine may result in proliferation with resultant changes in results with the protein test area. Urine from organisms metabolize the glucose. Contamination of the urine specimen may affect protein (and to a lesser extent

MATERIALS

- Material Provided**
- Strips
 - Pack


- Material required but not Provided**
- Specimen collection container

DIRECTIONS FOR USE

Allow the strip, urine specimen, (15-30°C) prior to testing.

1. Remove the strip from the closed canister. Immediately close the canister tightly. Completely immerse the reagent strip in the urine and immediately remove the strip to a clean surface below.
2. While removing the strip from the urine container to remove excess urine, bring the edge of the strip into contact with a dry towel to avoid mixing chemicals with urine. See illustration 2 below.
3. Compare the reagent areas to the color chart at the specified times. Hold the strip against the color chart. See illustration 3 below.

Note: Results may be read up to 2 minutes.

INTERPRETATIONS OF RESULTS

Results are obtained by direct comparison of the color blocks on the test strip with the color chart on the canister label. The color blocks represent nominal values. In the event of unexpected results, repeat the test using a new strip. Results are recommended; confirm that the date printed on the canister label, controls and repeat the test using a new strip immediately and contact your

QUALITY CONTROL

For best results, performance of reagent strips on positive and negative specimens/controls should be checked whenever a new canister is first opened and after the strip is opened for adequate standards of performance.

LIMITATIONS

Note: The Urinalysis Reagent Strips may give abnormal urine color such as drugs: Gantrisin®, Azo Gantano®, nitrofurantoin, riboflavin, etc. The color development of the test may be produced that could be interpreted as a positive result. **Ascorbic acid:** No interference is known.

Glucose: This test is highly specific for glucose. No substance excreted in urine other than glucose is known to give a positive result. The reagent area does not react with ketones, lactose, galactose, fructose or other metabolic substances, nor with reducing metabolites of drugs (e.g. salicylates and nalidixic acid). Sensitivity may be decreased in specimens with high specific gravity (>1.025) and with ascorbic acid concentrations of ≥ 25 mg/dL. High ketone levels ≥ 100 mg/dL may cause false negative results for specimens containing a small amount of glucose (50-100 mg/dL).

Bilirubin: Bilirubin is absent in normal urine, so any positive result, including a trace positive, indicates an underlying pathological condition and requires further investigation. Reactions may occur with urine containing large doses of chlorpromazine or rifampin that might be mistaken for positive bilirubin. The presence of bilirubin-derived bile pigments may mask the bilirubin reaction. This phenomenon is characterized by color development on the test patch that does not correlate with the colors on the color chart. Large concentrations of ascorbic acid may decrease sensitivity.

Ketone: The test does not react with acetone or β-hydroxybutyrate. Urine specimens of high pigment and other substances containing sulfhydryl groups occasionally give reactions up to and including trace (+).

higher more than 300 mg/dL may cause false positive results. Interference may be caused by non-ionic urine components such as urea, creatinine, etc. For accurate results, add 0.005 to the specific gravity reading. The presence of myoglobin, hemoglobin or intact erythrocytes may cause false positive results. Blue spots indicate intact erythrocytes. False positive results are provided for hemoglobin and for urea. False positive results are often seen with urine from menstruating women. High pH reduces sensitivity, while moderate pH reduces specificity. Microbial peroxidase, which is present in urine, may cause a false positive reaction. The test is not intended to detect intact erythrocytes.

Excess urine remains on the strip, which may cause a false positive result. Excess urine remains on the strip, in which the acid buffer from the protein test area may cause a false positive result. Excess urine remains on the strip, in which the acid buffer from the protein test area may cause a false positive result. Excess urine remains on the strip, in which the acid buffer from the protein test area may cause a false positive result.

Urobilinogen should be interpreted as absent if the result is negative. The absence of urobilinogen is not a reliable indicator of the absence of urobilinogen. Substances known to react with Ehrlich's reagent are urobilinogen, urobilin, and urobilinamides. False negative results may be obtained if the urine specimen is not used to detect porphobilinogen.

The test will not react with any other substance that may cause a false positive result. Form pink to red color should be interpreted as a positive result. Color intensity is not proportional to the amount of protein in the urine specimen. Pink spots or pink edges on the strip may indicate a positive result. Comparing the reacted reagent area on a color chart to the color chart may cause false negatives in urine containing urobilinogen. The sensitivity of this test is reduced for urine with a pH of 7 or higher.

For accurate results, antibiotics should be discontinued prior to testing. A negative result does not indicate the absence of urobilinogen. Negative results may occur in urinary porphobilinogen reductase to convert nitrite to nitrate; under a sufficient length of time (at least 4 hours) or when dietary nitrate is absent. Between 60-120 seconds to allow for complete reaction. The result is proportional to the number of protein molecules. High specific gravity or elevated glucose may cause false negative results. The test results to be artificially low. The high concentrations of oxalic acid may also cause false negative results. Penicillin may cause decreased reactivity, resulting in a false negative reaction. High urinary protein may cause false positive results. This test will not react with erythrocytes.

Journal of Clinical Science. CRC Crit. Rev. Clin. Chem. 1965.
Screening for Urinary Occult Blood, Protein, and Glucose in Urine. JAMA 201:129-132, 1967.
Regulation of Ketogenesis. J. Biol. Chem. 1965.

Why Ketone Bodies? Postgrad. Med. J. (June 1965) 41:372-378.
Concentrations in Plasma and Urine. Lancet. 1965.
Nitroprusside Test for Ketone Bodies in Urine. J. Clin. Chem. 1965.
Management by Laboratory Methods, 18th Ed. W.B. Saunders Company. 1976.
Practical Clinical Chemistry, 2nd Ed. 2205, 1994.

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BIBLIOGRAFIA

- Free AH, Free HM. *Urinalysis, Critical Discipline of Clinical Science*. CRC Crit. Rev. Clin. Lab. Sci. 3(4): 481-531, 1972.
- Yoder J, Adams EC, Free, AH. *Simultaneous Screening for Urinary Occult Blood, Protein, Glucose, and pH*. Amer. J. Med Tech. 31:285, 1965.
- Shcherbin B, Fritz H. *Subnormal Levels of Glucose in Urine*. JAMA 201:129-132, 1967.
- McGarry JD, Lilly, Lecture, 1978: New Perspectives in the Regulation of Ketogenesis. Diabetes 28: 517-523 May, 1978.
- Williamson DH. *Physiological Ketoses, or Why Ketone Bodies?* Postgrad. Med. J. (June Suppl.): 372-375, 1971.
- Paterson P et al. *Maternal and Fetal Ketone Concentrations in Plasma and Urine*. Lancet: 862-865; April 22, 1967.
- Fraser J, et al. *Studies with a Simplified Nitroprusside Test for Ketone Bodies in Urine, Serum, Plasma and Milk*. Clin. Chem. Acta II: 372-378, 1965.
- Henry JB, et al. *Clinical Diagnosis and Management by Laboratory Methods*, 18th Ed. Philadelphia. Saunders. 396-397, 415, 1991.
- Tietz NW. *Clinical Guide to Laboratory Tests*. W.B. Saunders Company. 1976.
- Burtis CA, Ashwood ER. *Tietz Textbook of Clinical Chemistry* 2nd Ed. 2205, 1994.

PRESENTACION

	Ref.	Cont	100 Tiras
52001-52002-52003-52005-52010-52011			