


Flexible Contact Slide for microbial screening of surfaces (ISO 18593)



Description	Packaging	Ref.
Contact Slide 1 (PCA + TTC + Neutralizing / V.R.B.G. Agar + Neutralizing) Total bacterial count, enterobacteria, with inactivation of disinfectants.	20 slides 120 slides	525262 53526
Contact Slide 2 (PCA + TTC + Neutralizing / Rose Bengal CAF Agar + Neutralizing) Total bacterial count, yeasts and moulds, with inactivation of disinfectants.	20 slides 120 slides	525272 53527
Contact Slide 3 (PCA+TTC + Neutralizing / Vogel Johnson Agar) Total bacterial count with inactivation of disinfectants, staphylococci.	20 slides 120 slides	525302 53530
Contact Slide 4 (Cetrimide Agar + Neutralizing / Rose Bengal CAF Agar + Neutralizing) <i>Pseudomonas</i> spp., yeasts and moulds.	20 slides 120 slides	525282 # 53528 #
Contact Slide 5 (V.R.B.G. Agar + Neutralizing / Bile Aesculin Agar + Neutralizing) Enterobacteria, fecal streptococci, with inactivation of disinfectants.	20 slides 120 slides	525312 53531
Contact Slide 8 (PCA + TTC + Neutralizing / V.R.B.L. Agar + Neutralizing) Total bacterial count, coliforms, with inactivation of disinfectants.	20 slides 120 slides	525622 53562
Contact Slide 9 (PCA + TTC + Neutralizing / PCA + TTC + Neutralizing) Total bacterial count with inactivation of disinfectants.	20 slides 120 slides	525322 53532
Contact Slide 10 (Rose Bengal CAF Agar + Neutralizing / Vogel Johnson Agar) Yeasts and moulds with inactivation of disinfectants, staphylococci.	20 slides 120 slides	525452 # 53545 #
Contact Slide Chrom 1 (Chromatic™ <i>E.coli</i> -coliform / Chromatic™ <i>Salmonella</i>) Chromogenic selective media for detecting <i>E. coli</i> , coliforms and <i>Salmonella</i> spp. on surfaces or in food and waters.	20 slides 120 slides	525292 53529
Contact Slide Chrom 2 (Chromatic™ <i>E.coli</i> -coliform / PCA+TTC + Neutralizing) Chromogenic selective medium for detecting <i>E. coli</i> , coliforms and for total bacterial count.	20 slides 120 slides	525392 53539
Contact Slide Chrom 3 (Chromatic™ <i>E.coli</i> -coliform / Baird Parker Agar + Neutralizing) Chromogenic selective medium for detecting <i>E. coli</i> , coliforms and <i>Staphylococcus aureus</i> with inactivation of disinfectants.	20 slides 120 slides	525462 # 53546 #
Contact Slide TSA + Neutralizing / TSA + Neutralizing Total bacterial count with inactivation of disinfectants.	20 slides 120 slides	525472 # 53547 #
Contact Slide TSA + Neutralizing / Rose Bengal CAF Agar + Neutralizing Total bacterial count with inactivation of disinfectants, yeasts and moulds.	20 slides 120 slides	525482 # 53548 #

Product with minimum order