

## BOÎTES DE PÉTRI / *PETRI DISHES*

DIAM. 90 mm - H 14.2 mm



Flèches imprimées sur sachet recyclable pour indiquer le sens d'ouverture / *Arrows printed on recyclable sleeve to indicate the opening way*

- ▶ 0/1/3 Ergots de ventilation
- ▶ Polystyrène cristal
- ▶ Planéité optimale
- ▶ Production classe 100

- ▶ 0/1/3 Vents
- ▶ Crystal Polystyrene
- ▶ Optimal flatness
- ▶ Production class 100


### Boîtes de Pétri / Petri dishes (diam. 90 mm H 14.2 mm)

Cat. No.	Ergots Vents	Pcs/gaine Units/sleeve	Gaines Sleeves	Sursac Outer bag	Qté/carton Qty/case	Sterile R SAL 10 <sup>-6</sup>
BP90B-03	0	33	25	2	825	-
BP90B-04	0	33	25	2	825	✓
SB90-01	0	33	25	0	825	-
BP91B-03	1	33	25	2	825	-
BP91B-04	1	33	25	2	825	✓
SB91-01	1	33	25	0	825	-
BP93B-15	3	33	25	2	825	-
BP93B-16	3	33	25	2	825	✓
BP93D-01	3	33	25	2	825	-
BP93D-02	3	33	25	2	825	✓
SB93-01	3	33	25	0	825	-
SB93-03	3	33	25	0	825	✓

  Zone de dépoli / Frosted area
  Sterile R SAL 10<sup>-6</sup>
 Aseptique / Aseptic

## 2 COMPARTIMENTS AVEC DÉPOLI / 2 COMPARTMENTS

DIAM. 90 mm - H 14.2 mm



- ▶ 3 ergots de ventilation
- ▶ Polystyrène cristal
- ▶ Planéité optimale
- ▶ Production classe 100
- ▶ 2 compartiments
- ▶ Zone de dépoli extérieure

- ▶ 3 vents
- ▶ Crystal Polystyrene
- ▶ Optimal flatness
- ▶ Production class100
- ▶ 2 compartiments
- ▶ Frosted area

### Boîtes de Pétri 2 compartiments / 2 compartment Petri dishes (diam. 90 mm H 14.2 mm)

Cat. No.	Ergots Vents	Pcs/gaine Units/sleeve	Gaines Sleeves	Sursac Outer bag	Qté/carton Qty/case	Sterile R SAL 10 <sup>-6</sup>
BI93D-01	3	33	25	2	825	-
BI93D-02	3	33	25	2	825	✓

  Sterile R SAL 10<sup>-6</sup>
 Aseptique / Aseptic
Indexage de la boîte pour couler les milieux dans leur compartiment respectif / *Allows indexing to facilitate automated pouring of media into their respective compartments*