

PIPETTES PAILLE / STRAW PIPETS

1 mL -> 5 mL




- ▶ Pipettes paille en Polypropylène
- ▶ Conditionnées en sachet ZIP avec témoin d'intégrité
- ▶ Prélèvement de liquide
- ▶ Dimension universelle
- ▶ Traçabilité sur chaque sachet
- ▶ Straw pipets in Polypropylene
- ▶ Packed in tamper-evident ZIP bags
- ▶ Liquid sampling
- ▶ Universal size
- ▶ Traceability on each bag

Pipettes paille en PP / Straw pipets in PP


Cat. No.	Volume (mL)	L (mm)	Emballage	Qté/boîte Qty/box	Qté/carton Qty/case	Sterile SAL 10 ⁻³
EPI1B-01	1	190	25	40x12	12000	✓
EPI1B-05	1	150	100	10x12	12000	✓
EPI24B-01	2.2	240	25	240	6000	✓
EPI24B-02	2.2	240	50	120	6000	✓
EPI5-01	5	187	10	30x12	3600	✓

PIPETTES DE TRANSFERT / TRANSFER PIPETS

1 mL & 3 mL



- ▶ Pipettes graduées en Polyéthylène
- ▶ Souple
- ▶ Ecoulement rapide et total
- ▶ Traçabilité sur chaque sachet
- ▶ Graduated pipets in Polyethylene
- ▶ Supple
- ▶ Complete and quick removal of liquid
- ▶ Traceability on each bag



Pipettes de transfert en PE / Transfer pipets in PE

Cat. No.	Volume (mL)	L (mm)	Emballage	Qté/boîte Qty/box	Qté/carton Qty/case	Sterile SAL 10 ⁻³
PAS1-00	1	150	-	500x10	5000	-
PAS1B-01	1	225	20	200x10	2000	✓
PAS1A-01	1	150	1	200x10	2000	✓
PAS3-00	3	150	-	500x10	5000	-
PAS3A-01	3	150	1	200x10	2000	✓

Sachet pelable
Peelable bag