

5

DEDICATED LINES  
*LÍNEAS DEDICADAS*



## SPINLAB / SELECTRA - FLEXOR SYSTEM

PRODUCT · PRODUCTO	METHOD · MÉTODO	SIZE · PRESENTACIÓN	REF. · CÓDIGO
BIOCHEMISTRY · BIOQUÍMICA			
ALBUMIN · ALBÚMINA	Bromocresol green. Colorimetric	10x25 mL	SP1001020
URINE & CSF TOTAL PROTEINS · PROTEÍNAS TOTALES EN ORINA Y LCR	Red Pyrogallol. Colorimetric	10x25 mL	SP1001025
TOTAL BILIRUBIN · BILIRRUBINA TOTAL	DMSO. Colorimetric	10x25 / 2x5 mL	SP1001042
DIRECT BILIRUBIN · BILIRUBINA DIRECTA	DMSO. Colorimetric	10x25 / 2x5 mL	SP1001043
BILIRUBIN T&D · BILIRUBINA TyD	DMSO. Colorimetric	10x25 / 2x5 mL	SP1001044
CALCIUM · CALCIO	v/v o-Cresolphthalein. Colorimetric	5x25 / 5x25 mL	SP1001061
CALCIUM · CALCIO	Arsenazo III Colorimetric	10x25 mL	SP1001065
HDL DIRECT · DIRECTO	Direct. Colorimetric	10x24 mL / 10x8 mL	SP1001096
HDL DIRECT · DIRECTO	Direct. Colorimetric	2x24 mL / 2x8 mL	SP1001097
CREATININE · CREATININA	Jaffé Colorimetric-Kinetic	5x25 / 5x25 mL	SP1001111
CREATININE · CREATININA	Trinder Colorimetric	10x24 mL / 10x8 mL	SP1001117
PHOSPHORUS · FÓSFORO	Phosphomolybdate. UV	10x25 mL	SP1001155
IRON · HIERRO	Ferrozine Colorimetric	10x25 mL / 3x5 mL	SP1001247
LIPASE · LIPASA	Enzymatic. Kinetic - Colorimetric	2x20 mL / 2x4 mL 10x20 mL / 5x8 mL	SP1001274
MAGNESIUM · MAGNESIO	Xylidyl Blue. Colorimetric	10x25 mL	SP1001285
TOTAL PROTEINS · PROTEÍNAS TOTALES	Biuret. Colorimetric	10x25 mL	SP1001291
CHLORIDE · CLORURO	Thiocyanate-Hg. Color	10x25 mL	SP1001360
URIC ACID · ÁCIDO ÚRICO	Uricase-PAP. Enz-Color	5x25 mL / 5x25 mL	SP41001
GLUCOSE · GLUCOSA	GOD-PAP Enzymatic Colorimetric	10x25 mL	SP41011
CHOLESTEROL · COLESTEROL	CHOD/PAP Enzymatic Colorimetric	10x25 mL	SP41021
LDL DIRECT · DIRECTO	Enzymatic Colorimetric	10x24 mL / 10x8 mL 2x24 mL / 2x8 mL	SP41023 SP41024
TRIGLYCERIDES · TRIGLICÉRIDOS	GPO-POD. Enzymatic Colorimetric	10x25 mL	SP41031
UREA	Urease - GLDH. Kinetic	10x20 mL / 10x5 mL	SP41041
AMYLASE · AMILASA	CNPG3 Kinetic	10x25 mL	SP41201
LDH	Kinetic - UV	10x20 mL / 10x5 mL	SP41214
ALP	DGKC. Kinetic optimized	10x20 mL / 10x5 mL	SP41233
CK-NAC	NAC Immunoinh. Kinetic-UV	10x20 mL / 10x5 mL	SP41250
CK-MB	Anti CK-M Immunoinh. Kin.-UV	10x20 mL / 10x5 mL	SP41254
GOT / AST	IFCC Enzymatic - UV	10x20 mL / 10x5 mL	SP41264
GPT / ALT	IFCC Enzymatic - UV	10x20 mL / 10x5 mL	SP41274
γ-GT	Carboxy substrate Kinetic	10x20 mL / 10x5 mL	SP41288



PRODUCT · PRODUCTO	METHOD · MÉTODO	SIZE · PRESENTACIÓN	REF. · CÓDIGO
<b>IMMUNOTURBIDIMETRY / TURBILATEX</b>			
<b>APOLIPOPROTEIN A-I</b>	Immunoturbidimetry	2x20 mL / 2x5 mL	SP1103012
<b>APOLIPOPROTEIN B</b>	Immunoturbidimetry	2x20 mL / 2x5 mL	SP1103013
<b>ANTITHROMBIN III</b>	Immunoturbidimetry	2x20 mL / 2x5 mL	SP1102014
<b>α-1-ANTITRYPSIN</b>	Immunoturbidimetry	2x20 mL / 2x5 mL	SP1102054
<b>TRANSFERRIN</b>	Immunoturbidimetry	2x20 mL / 2x5 mL	SP1102134
<b>α-1-GLYCOPROTEIN ACID</b>	Immunoturbidimetry	2x20 mL / 2x5 mL	SP1102154
<b>IMMUNOGLOBULIN G (IgG)</b>	Immunoturbidimetry	2x20 mL / 2x5 mL	SP1103004
<b>IMMUNOGLOBULIN A (IgA)</b>	Immunoturbidimetry	2x20 mL / 2x5 mL	SP1103014
<b>IMMUNOGLOBULIN M (IgM)</b>	Immunoturbidimetry	2x20 mL / 2x5 mL	SP1103024
<b>CRP</b>	Turbilatex	2x20 mL / 2x5 mL	SP1107001
<b>RF</b>	Turbilatex	2x20 mL / 2x5 mL	SP1107005
<b>ASO</b>	Turbilatex	2x20 mL / 2x5 mL	SP1107010
<b>Lp(a)</b>	Turbilatex	1x20 mL / 1x4 mL	SP1107020
<b>β-2-MICROGLOBULIN</b>	Turbilatex	2x20mL / 2x5 mL	SP1107030
<b>FERRITIN</b>	Turbilatex	2x20 mL / 2x5 mL	SP1107040
<b>IMMUNOGLOBULIN E (IgE)</b>	Turbilatex	2x10 mL / 2x5 mL	SP1107050
<b>MICROALBUMIN</b>	Turbilatex	2x20 mL / 2x5 mL	SP1107070

## TOKYO BOEKI / SPINTECH 240 / BIOLIS 24I / PREMIUM 24I SYSTEM

PRODUCT · PRODUCTO	METHOD · MÉTODO	SIZE · PRESENTACIÓN	REF. · CÓDIGO
BIOCHEMISTRY · BIOQUÍMICA			
ALBUMIN · ALBÚMINA	Bromocresol green colorimetric	10x35 mL	TK1001020
URINE & CSF TOTAL PROTEINS · PROTEÍNAS TOTALES EN ORINA Y LCR	Red Pyrogallol. Colorimetric	10x35 mL	TK1001025
TOTAL BILIRUBIN · BILIRRUBINA TOTAL	DMSO. Colorimetric	8x35 mL / 1x10 mL	TK1001042
DIRECT BILIRUBIN · BILIRUBINA DIRECTA	DMSO. Colorimetric	8x35 mL / 1x10 mL	TK1001043
CALCIUM · CALCIO	Arsenazo III Colorimetric	10x35 mL	TK1001065
HDL DIRECT · DIRECTO	Direct. Colorimetric	10x24 mL / 10x8 mL	TK1001096
CREATININE · CREATININA	Jaffé Colorimetric-Kinetic	8x20 mL / 12x13 mL	TK1001111
CREATININE · CREATININA	Trinder Colorimetric	10x24 mL / 10x8 mL	TK1001117
CHOLINESTERASE	Butyrylthiocholine	2x25 mL / 2x5 mL	TK41210
PHOSPHORUS · FÓSFORO	Phosphomolybdate. UV	10x35 mL	TK1001155
IRON · HIERRO	Ferrozine Colorimetric	9x24 mL / 2x10 mL	TK1001247
LIPASE · LIPASA	Enzymatic. Kinetic - Colorimetric	2x25 mL / 2x5 mL	TK1001275
MAGNESIUM · MAGNESIO	Xylidyl Blue. Colorimetric	10x35 mL	TK1001285
TOTAL PROTEINS · PROTEÍNAS TOTALES	Biuret. Colorimetric	10x35 mL	TK1001291
CHLORIDE · CLORURO	Thiocyanate-Hg. Color	10x35 mL	TK1001360
URIC ACID · ÁCIDO ÚRICO	Uricase-PAP. Enz-Color	8x20 mL / 12x13 mL	TK41001
GLUCOSE · GLUCOSA	GOD-PAP Enzymatic Colorimetric	10x35 mL	TK41011
CHOLESTEROL · COLESTEROL	CHOD/PAP Enzymatic Color.	10x35 mL	TK41021
LDL DIRECT · DIRECTO	Enzymatic Colorimetric	10x24 mL / 10x8 mL	TK41023
TRIGLYCERIDES · TRIGLICÉRIDOS	GPO-POD. Enz.-Colorimetric	10x35 mL	TK41031
UREA	Urease - GLDH. Kinetic	10x25 mL / 10x7 mL	TK41041
AMYLASE · AMILASA	CNP3. Kinetic	10x35 mL	TK41201
		2x35 mL	TK41202
LDH	IFCC. Kinetic - UV	10x25 mL / 10x7 mL	TK41214
ALP	DGKC. Kinetic optimized	10x25 mL / 10x7 mL	TK41233
CK-NAC	NAC Immunoinh. Kinetic-UV	10x25 mL / 10x7 mL	TK41250
CK-MB	Anti CK-M Immunoinh. Kinetic-UV	10x25 mL / 10x7 mL	TK41254
GOT / AST	IFCC Enzymatic - UV	10x25 mL / 10x7 mL	TK41264
GPT / ALT	IFCC Enzymatic - UV	10x25 mL / 10x7 mL	TK41274
γ-GT	Carboxy substrate Kinetic	10x25 mL / 10x7 mL	TK41288



PRODUCT · PRODUCTO	METHOD · MÉTODO	SIZE · PRESENTACIÓN	REF. · CÓDIGO
IMMUNOTURBIDIMETRY / TURBILATEX			
APOLIPOPROTEIN A-I	Immunoturbidimetry	2x24 mL / 2x6 mL	TK1103012
APOLIPOPROTEIN B	Immunoturbidimetry	2x24 mL / 2x6 mL	TK1103013
ANTITHROMBIN III	Immunoturbidimetry	2x24 mL / 2x6 mL	TK1102014
α-1-ANTITRYPSIN	Immunoturbidimetry	2x24 mL / 2x6 mL	TK1102054
CERULOPLASMIN	Immunoturbidimetry	2x24 mL / 2x6 mL	TK1102064
HAPTOGLOBIN	Immunoturbidimetry	2x24 mL / 2x6 mL	TK1102074
COMPLEMENT C3	Immunoturbidimetry	2x24 mL / 2x6 mL	TK1102094
COMPLEMENT C4	Immunoturbidimetry	2x24 mL / 2x6 mL	TK1102104
PREALBUMIN	Immunoturbidimetry	2x24 mL / 2x6 mL	TK1102124
TRANSFERRIN	Immunoturbidimetry	2x24 mL / 2x6 mL	TK1102134
α-1-GLYCOPROTEIN ACID	Immunoturbidimetry	2x24 mL / 2x6 mL	TK1102154
IMMUNOGLOBULIN G (IgG)	Immunoturbidimetry	2x24 mL / 2x6 mL	TK1103004
IMMUNOGLOBULIN A (IgA)	Immunoturbidimetry	2x24 mL / 2x6 mL	TK1103014
IMMUNOGLOBULIN M (IgM)	Immunoturbid	2x24 mL / 2x6 mL	TK1103024
CRP	Turbilatex	2x24 mL / 1x12 mL	TK1107001
RF	Turbilatex	2x24 mL / 1x12 mL	TK1107005
ASO	Turbilatex	2x24 mL / 1x12 mL	TK1107010
Lp(a)	Turbilatex	2x10 mL / 1x4 mL	TK1107020
β-2-MICROGLOBULIN	Turbilatex	2x24 mL / 1x12 mL / 1x1 mL (Cal)	TK1107030
FERRITIN	Turbilatex	1x18 mL / 1x10 mL	TK1107040
		2x24 mL / 1x12 mL	TK1107140
IMMUNOGLOBULIN E (IgE)	Turbilatex	2x10 mL / 2x5 mL	TK1107050
MYOGLOBIN	Turbilatex	2x10 mL / 2x5 mL	TK1107060
MICROALBUMIN	Turbilatex	2x24 mL / 1x12 mL / 1x1 mL (Cal)	TK1107070
Glycated Haemoglobin (HbA1c)	Turbilatex	1x25 mL (200 Tests)	TK43090

## SPIN 120 / MINDRAY BS 120 / BS 200 SYSTEM

PRODUCT · PRODUCTO	METHOD · MÉTODO	SIZE · PRESENTACIÓN	REF. · CÓDIGO
BIOCHEMISTRY · BIOQUÍMICA			
ALBUMIN · ALBÚMINA	Bromocresol green. Colorimetric	6x30 mL	MI1001020
URINE & CSF TOTAL PROTEINS · PROTEÍNAS TOTALES EN ORINA Y LCR	Red Pyrogallol. Colorimetric	6x30 mL	MI1001025
TOTAL BILIRUBIN · BILIRRUBINA TOTAL	DMSO. Colorimetric	5x30 mL / 1x10 mL	MI1001042
DIRECT BILIRUBIN · BILIRUBINA DIRECTA	DMSO. Colorimetric	5x30 mL / 1x10 mL	MI1001043
CALCIUM · CALCIO	Arsenazo III Colorimetric	6x30 mL	MI1001065
HDL DIRECT · DIRECTO	Direct. Colorimetric	4x30 mL / 2x20 mL	MI1001096
CREATININE · CREATININA	Jaffé Colorimetric-Kinetic	3x30 mL / 3x30 mL	MI1001111
CREATININE · CREATININA	Trinder Colorimetric	4x30 mL / 2x20 mL	MI1001117
PHOSPHORUS · FÓSFORO	Phosphomolybdate. UV	6x30 mL	MI1001155
IRON · HIERRO	Ferrozine Colorimetric	4x30 mL / 1x10 mL	MI1001247
LIPASE · LIPASA	Enzymatic. Kinetic - Colorimetric	2x30 mL / 1x15 mL	MI1001275
MAGNESIUM · MAGNESIO	Xylidyl Blue. Colorimetric	6x30 mL	MI1001285
TOTAL PROTEINS · PROTEÍNAS TOTALES	Biuret. Colorimetric	6x30 mL	MI1001291
CHLORIDE · CLORURO	Thiocyanate-Hg. Colorimetric	6x30 mL	MI1001360
URIC ACID · ÁCIDO ÚRICO	Uricase-PAP. Enz-Colorimetric	3x30 mL / 3x30 mL	MI41001
GLUCOSE · GLUCOSA	GOD-PAP Enzymatic Colorimetric	6x30 mL	MI41011
CHOLESTEROL · COLESTEROL	CHOD/PAP Enzymatic Colorimetric	6x30 mL	MI41021
LDL DIRECT · DIRECTO	Enzymatic Colorimetric	4x30 mL / 2x20 mL	MI41023
TRIGLYCERIDES · TRIGLICÉRIDOS	GPO-POD. Enz.-Colorimetric	6x30 mL	MI41031
UREA	Urease - GLDH. Kinetic	5x25 mL / 1x32 mL	MI41041
AMYLASE · AMILASA	CNPG3. Kinetic	6x30 mL	MI41201
LDH	IFCC. Kinetic - UV	5x25 mL / 1x32 mL	MI41214
ALP	DGKC. Kinetic optimized	5x25 mL / 1x32 mL	MI41233
GOT / AST	IFCC Enzymatic - UV	5x25 mL / 1x32 mL	MI41264
GPT / ALT	IFCC Enzymatic - UV	5x25 mL / 1x32 mL	MI41274
γ-GT	Carboxy substrate Kinetic	5x25 mL / 1x32 mL	MI41288



PRODUCT · PRODUCTO	METHOD · MÉTODO	SIZE · PRESENTACIÓN	REF. · CÓDIGO
<b>IMMUNOTURBIDIMETRY / TURBILATEX</b>			
<b>APOLIPOPROTEIN A-I</b>	Immunoturbidimetry	2x30 mL / 1x15 mL	MI1103012
<b>APOLIPOPROTEIN B</b>	Immunoturbidimetry	2x30 mL / 1x15 mL	MI1103013
<b>ANTITHROMBIN III</b>	Immunoturbidimetry	2x30 mL / 1x15 mL	MI1102014
<b>α-1-ANTITRYPSIN</b>	Immunoturbidimetry	2x30 mL / 1x15 mL	MI1102054
<b>CERULOPLASMIN</b>	Immunoturbidimetry	2x30 mL / 1x15 mL	MI1102064
<b>HAPTOGLOBIN</b>	Immunoturbidimetry	2x30 mL / 1x15 mL	MI1102074
<b>COMPLEMENT C3</b>	Immunoturbidimetry	2x30 mL / 1x15 mL	MI1102094
<b>COMPLEMENT C4</b>	Immunoturbidimetry	2x30 mL / 1x15 mL	MI1102104
<b>TRANSFERRIN</b>	Immunoturbidimetry	2x30 mL / 1x15 mL	MI1102134
<b>α-1-GLYCOPROTEIN ACID</b>	Immunoturbidimetry	2x30 mL / 1x15 mL	MI1102154
<b>IMMUNOGLOBULIN G (IgG)</b>	Immunoturbidimetry	2x30 mL / 1x15 mL	MI1103004
<b>IMMUNOGLOBULIN A (IgA)</b>	Immunoturbidimetry	2x30 mL / 1x15 mL	MI1103014
<b>IMMUNOGLOBULIN M (IgM)</b>	Immunoturbidimetry	2x30 mL / 1x15 mL	MI1103024
<b>CRP</b>	Turbilatex	2x30 mL / 1x15 mL	MI1107001
<b>RF</b>	Turbilatex	2x30 mL / 1x15 mL	MI1107005
<b>ASO</b>	Turbilatex	2x30 mL / 1x15 mL	MI1107010
<b>Lp(a)</b>	Turbilatex	1x20 mL / 1x4 mL	MI1107020
<b>β-2-MICROGLOBULIN</b>	Turbilatex	2x30 mL / 1x15 mL	MI1107030
<b>FERRITIN</b>	Turbilatex	2x30 mL / 1x15 mL	MI1107040
<b>IMMUNOGLOBULIN E (IgE)</b>	Turbilatex	2x10 mL / 2x5 mL	MI1107050
<b>MICROALBUMIN</b>	Turbilatex	2x30 mL / 1x15 mL	MI1107070



## MINDRAY BS 300 / BS 380 / BS 400 / SPIN 640 SYSTEM

PRODUCT · PRODUCTO	METHOD · MÉTODO	SIZE · PRESENTACIÓN	REF. · CÓDIGO
BIOCHEMISTRY · BIOQUÍMICA			
ALBUMIN · ALBÚMINA	Bromocresol green colorimetric	6x40 mL	MD1001020
URINE & CSF TOTAL PROTEINS · PROTEÍNAS TOTALES EN ORINA Y LCR	Red Pyrogallol. Colorimetric	6x40 mL	MD1001025
TOTAL BILIRUBIN · BILIRRUBINA TOTAL	DMSO. Colorimetric	5x40 mL / 1x10 mL	MD1001042
DIRECT BILIRUBIN · BILIRUBINA DIRECTA	DMSO. Colorimetric	5x40 mL / 1x10 mL	MD1001043
CALCIUM · CALCIO	Arsenazo III Colorimetric	6x40 mL	MD1001065
HDL DIRECT · DIRECTO	Direct. Colorimetric	4x30 mL / 2x20 mL	MD1001096
CREATININE · CREATININA	Jaffé Colorimetric-Kinetic	3x40 mL / 3x40 mL	MD1001111
CREATININE · CREATININA	Trinder Colorimetric	4x30 mL / 2x20 mL	MD1001117
PHOSPHORUS · FÓSFORO	Phosphomolybdate. UV	6x40 mL	MD1001155
IRON · HIERRO	Ferrozine Colorimetric	4x30 mL / 1x10 mL	MD1001247
LIPASE · LIPASA	Enzymatic. Kinetic - Colorimetric	2x40 mL / 2x8 mL	MD1001275
MAGNESIUM · MAGNESIO	Xylidyl Blue. Colorimetric	6x40 mL	MD1001285
TOTAL PROTEINS · PROTEÍNAS TOTALES	Biuret. Colorimetric	6x40 mL	MD1001291
CHLORIDE · CLORURO	Thiocyanate-Hg. Color	6x40 mL	MD1001360
URIC ACID · ÁCIDO ÚRICO	Uricase-PAP. Enz-Color	3x40 mL / 3x40 mL	MD41001
GLUCOSE · GLUCOSA	GOD-PAP Enzymatic Colorimetric	6x40 mL	MD41011
CHOLESTEROL · COLESTEROL	CHOD/PAP Enzymatic Color.	6x40 mL	MD41021
LDL DIRECT · DIRECTO	Enzymatic Colorimetric	4x30 mL / 2x20 mL	MD41023
TRIGLYCERIDES · TRIGLICÉRIDOS	GPO-POD. Enz.-Colorimetric	6x40 mL	MD41031
UREA	Urease-GLDH. Kinetic	4x40 mL / 2x20 mL	MD41041
AMYLASE · AMILASA	CNPG3. Kinetic	6x40 mL	MD41201
LDH	IFCC. Kinetic - UV	4x40 mL / 2x20 mL	MD41214
ALP	DGKC. Kinetic optimized	4x40 mL / 2x20 mL	MD41233
CK-NAC	NAC Immunoinh. Kinetic-UV	4x40 mL / 2x20 mL	MD41250
CK-MB	Anti CK-M Immunoinh. Kin.-UV	4x40 mL / 2x20 mL	MD41254
GOT / AST	IFCC Enzymatic - UV	4x40 mL / 2x20 mL	MD41264
GPT / ALT	IFCC. Enzymatic - UV	4X40 mL / 2x20 mL	MD41274
γ-GT	Carboxy substrate Kinetic	4X40 mL / 2x20 mL	MD41288





PRODUCT · PRODUCTO	METHOD · MÉTODO	SIZE · PRESENTACIÓN	REF. · CÓDIGO
<b>IMMUNOTURBIDIMETRY / TURBILATEX</b>			
APOLIPOPROTEIN A-I	Immunoturbidimetry	2x40 mL / 1x20 mL	MD1103012
APOLIPOPROTEIN B	Immunoturbidimetry	2x40 mL / 1x20 mL	MD1103013
ANTITHROMBIN III	Immunoturbidimetry	2x40 mL / 1x20 mL	MD1102014
$\alpha$ -1-ANTITRYPSIN	Immunoturbidimetry	2x40 mL / 1x20 mL	MD1102054
CERULOPLASMIN	Immunoturbidimetry	2x40 mL / 1x20 mL	MD1102064
HAPTOGLOBIN	Immunoturbidimetry	2x40 mL / 1x20 mL	MD1102074
COMPLEMENT C3	Immunoturbidimetry	2x40 mL / 1x20 mL	MD1102094
COMPLEMENT C4	Immunoturbidimetry	2x40 mL / 1x20 mL	MD1102104
TRANSFERRIN	Immunoturbidimetry	2x40 mL / 1x20 mL	MD1102134
$\alpha$ -1-GLYCOPROTEIN ACID	Immunoturbidimetry	2x40 mL / 1x20 mL	MD1102154
IMMUNOGLOBULIN G (IgG)	Immunoturbidimetry	2x40 mL / 1x20 mL	MD1103004
IMMUNOGLOBULIN A (IgA)	Immunoturbidimetry	2x40 mL / 1x20 mL	MD1103014
IMMUNOGLOBULIN M (IgM)	Immunoturbidimetry	2x40 mL / 1x20 mL	MD1103024
CRP	Turbilatex	2x40 mL / 1x20 mL	MD1107101
RF	Turbilatex	2x40 mL / 1x20 mL	MD1107105
ASO	Turbilatex	2x40 mL / 1x20 mL	MD1107110
Lp(a)	Turbilatex	2x40 mL / 2x8 mL	MD1107020
$\beta$ -2-MICROGLOBULIN	Turbilatex	2x40 mL / 1x20 mL	MD1107130
FERRITIN	Turbilatex	2x40 mL / 1x20 mL	MD1107140
IMMUNOGLOBULIN E (IgE)	Turbilatex	2x40 mL / 2x20 mL	MD1107050
MICROALBUMIN	Turbilatex	2x40 mL / 1x20 mL	MD1107170
CRP ULTRA	Turbilatex	2x40 mL / 1x20 mL	MD43134
<b>DRUGS OF ABUSE · DROGAS DE ABUSO</b>			
AMPHETAMINE EIA Urine	Kinetic UV	2x40 mL / 2x15 mL	MD930008
BENZODIACEPINES EIA Urine	Kinetic UV	2x40 mL / 2x15 mL	MD930018
CANNABINOID EIA Urine	Kinetic UV	2x40 mL / 2x15 mL	MD930028
COCAINE Metabolite EIA Urine	Kinetic UV	2x40 mL / 2x15 mL	MD930038
COTININE EIA Urine	Kinetic UV	2x40 mL / 2x15 mL	MD930048
ECSTASY EIA Urine	Kinetic UV	2x40 mL / 2x15 mL	MD930058
METHADONE EIA Urine	Kinetic UV	2x40 mL / 2x15 mL	MD930068
EDDP EIA Urine	Kinetic UV	2x40 mL / 2x15 mL	MD930078
OPIATES EIA Urine	Kinetic UV	2x40 mL / 2x15 mL	MD930088
BARBITURATES EIA Urine	Kinetic UV	2x40 mL / 2x15 mL	MD930098
ETHYL ALCOHOL EIA Urine	Kinetic UV	2x40 mL / 2x15 mL	MD930108

