

Dip-Slide for the screening of Urinary Infections

Product with minimum order



Description	Packaging	Ref.
URITEST (Mac Conkey Agar / Cetrimide Agar / Cled Andrade Agar) Enterobacteria, <i>Pseudomonas</i> spp., total bacterial count.	20 slides 120 slides	500152 51015
URITEST N (Mac Conkey Agar / Cetrimide Agar / Cled Agar) Enterobacteria, <i>Pseudomonas</i> spp., total bacterial count.	20 slides 120 slides 500 slides	500232 51023 51123
URITEST C (Mac Conkey MUG Agar / Cetrimide Agar / Cled Andrade Agar) Enterobacteria, <i>Escherichia coli</i> , <i>Pseudomonas</i> spp., total bacterial count.	20 slides 120 slides	500242 # 51024 #
URITEST EC (Mac Conkey Agar / E.COLI Agar / Cled Agar) Enterobacteria, <i>Escherichia coli</i> and coliforms, total bacterial count.	20 slides 120 slides	500412 51041
URITEST 2 (Cled Agar / Mac Conkey Agar) Total bacterial count, enterobacteria.	20 slides 120 slides 500 slides	500302 51030 51130
URITEST M (Mac Conkey Agar / Malt SELECTIVE Agar / Cled Andrade Agar) Enterobacteria, yeasts, total bacterial count.	20 slides 120 slides 500 slides	500182 51018 51118
URITEST MALTO (Mac Conkey Agar / Malt Agar / Cled Agar) Enterobacteria, yeasts, total bacterial count.	20 slides 120 slides 500 slides	500402 51040 51140
URITEST EF (Mac Conkey Agar / Bile Aesculin Agar / Cled Agar) Enterobacteria, enterococci, total bacterial count.	20 slides 120 slides	500702 51070
URITEST PENTA (Cled Agar / Mac Conkey Agar / Cetrimide Agar / Malt Agar / Bile Aesculin Agar) Total bacterial count, enterobacteria, <i>Pseudomonas</i> spp., yeasts, enterococci.	20 slides 120 slides	500142 51014

Dip-Slide for the screening of Urogenital Infections



Description	Packaging	Ref.
VAGITEST (Chocolate / Rogosa / Thayer Martin) Fastidious microorganisms, lactobacilli, <i>Neisseria gonorrhoeae</i> .	10 slides 120 slides	50020 51020

Dip-Slide for the screening of Dermatological Infections

Description	Packaging	Ref.
DERMATEST (Malto / Corn Meal / DTM-Taplin) Pathogenic dermatophytes and others non pathogenic fungi.	10 slides 120 slides	50021 51021